



"High performance sieves by masse präzise"



+91 6264915873
+91 6264901140



sales@masseprazise.com
www.masseprazise.com




Test Sieve



Available sizes in compliance with ISO 3310-1:2016/ASTM E11-22/BAA 410:1986

Nominal Aperture Supplementary Size (R40/3)	2.36	2	1.7	1.18	1	850	710	600	500	425	355	250	212	180	150	125	106	90	75	63	45	38	
	mm	mm	mm	mm	mm	µm	µm	µm	µm	µm	µm	µm	µm	µm	µm	µm	µm	µm	µm	µm	µm	µm	µm
Equivalent ASTM Number	8	10	12	16	18	20	25	30	35	40	45	60	70	80	100	120	140	170	200	230	325	400	
Equivalent BSS Number	7	8	10	14	16	18	22	25	30	36	44	52	72	85	100	120	150	170	200	240	350	400	
Nominal Wire Diameter (± 15%)	1 mm	0.9 mm	0.8 mm	0.63 mm	0.56 mm	500 µm	450 µm	400 µm	315 µm	280 µm	224 µm	160 µm	140 µm	125 µm	100 µm	90 µm	71 µm	63 µm	50 µm	45 µm	32 µm	30 µm	

Note: Sieve sizes 4mm, 3.35mm, 2.80mm & 32 µm made available as per requirement.



pen microscope



Linear Mesh Counter



Sieves Shaker



nominal aperture calibrator



Unmatched Quality & Durability

Premium Test Sieves & sieve accessories for Every Application

Engineered to Meet Industry Standards and Specifications

Wide Range of Sizes and Mesh Specifications